


Using an in-situ spectrophotometer to study dissolved organic carbon in a New England forest



Harry Jin^{1,2}, Abigail Whittington², Cédric G. Fichot²

Middleton High School, 2100 Bristol St, Middleton, WI 53562¹,

Department of Earth & Environment, Boston University, 685 Commonwealth Avenue, Boston, MA 02215²

Introduction

- Harvard Forest Long-Term Ecological Reserve (HF) is a forest research site in central MA
- Forest **carbon cycling** is key to global carbon models and understanding the impacts of climate change
- Attention to stream water **dissolved organic carbon** (DOC) and **dissolved organic matter** (DOM) because streams transport carbon between carbon pools (plant matter, soil, etc.)

- Traditional methods of studying DOC via grab sampling have drawbacks
- **Colored dissolved organic matter** (CDOM) is the light-absorbing portion of DOC and is closely tied to overall DOC
- Measurements of CDOM facilitate study of DOC

PURPOSE OF THIS STUDY:

- Enable continuous in situ measurement to study CDOM/DOC with greater detail and accuracy
- Understand DOC dynamics in forested wetlands and assess future climate impacts on carbon cycle



Fig. 1 Example of CDOM in Harvard Forest stream water samples. The different colors are due to varying amounts of CDOM.

Study site and methods

- Studied a stream in HF named Arthur Brook Lower
- Drains Beaver Swamp, a recently formed swamp
- Beaver Swamp is fed by Arthur Brook Upper
- Used **UV-Visible spectrophotometry** to study CDOM
- Measures across range of wavelengths to produce CDOM absorption spectrum
- **Absorption coefficient a_g** quantifies the amount of light absorbed by CDOM in a sample



Fig. 3 s::can spectro::lyser set up in Arthur Lower.



Fig. 2 Laboratory spectrophotometer.

- Magnitude and shape of CDOM spectrum related to DOC concentration and composition
- Spectrophotometry can be done in the lab (Fig. 2) or in-situ (Fig 3.)
- Using s::can spectro::lyser for in-situ measurements
- Sensor deployed during fall 2023 and spring 2024

Correction scheme




Fig. 4 Using a correction scheme to account for errors improves CDOM spectra retrieval.

- An **empirical correction scheme** was developed and applied to in-situ spectra from the s::can
- Lab CDOM spectra and measurements used as reference
- Scheme corrects for errors due to sensor drift, detection limits, instrument features, etc.




Fig. 5 CDOM Absorption at 250 nm, laboratory measurements and corrected in-situ data. Avg % error = < 1%

DOC estimation

- CDOM absorption is closely tied to overall DOM/DOC properties and can be used as a **proxy for DOC concentration**




Fig. 6 Quadratic fit between DOC and CDOM absorption coefficient at 250 nm ($r^2 = 0.967$). Applied to all s::can data.

- Relationship between DOC and absorption differs between environments
- DOC **varies nonlinearly** with a_g in Arthur Brook
- Predicted with a **quadratic fit**
- Model using s::can absorption predicts DOC concentration within +/- 5.9%
- Error distributed evenly around 0




Fig. 7 Evaluation of DOC prediction. Mean absolute pct. error = 5.88%.

Calculation of DOC flux

- Continuous s::can measurements yield high resolution DOC data
- **DOC flux**, the mass of DOC moved by the stream per unit time, was calculated as: DOC concentration * stream discharge




Fig. 8 DOC and hydrological data for fall 2023 and spring 2024. Discharge and precipitation data from the HF data archive.

- Total fall DOC flux: 1021 kg; total spring DOC flux: 1017 kg
- Compared to 947 kg (fall) and 913 kg (spring) estimated from grab samples
- **Annual DOC export**: approx. 4 metric tons




Fig. 9 Cumulative DOC fluxes based on in-situ observations and grab samples. In-situ estimates of total DOC were 7.8% higher in fall, and 11.4% higher in spring.

Discussion

- Absorption spectra from the s::can spectro::lyser can be **successfully corrected** to match laboratory measurements
- Correction crucial for measuring advanced parameters such as spectral slope; useful for more detailed DOC analysis in the future
- **Future deployment** of s::can instrument to more remote environments is feasible
- Continuous observation and high temporal resolution are valuable
- Detection of short but extreme “pulses” of DOC improved, giving more accurate flux estimates
- **Discharge** is closely tied to DOC concentration and is the **main driver of DOC flux**
- Suggests Arthur Lower’s DOC is largely material flushed from Beaver Swamp
- Study of Arthur Upper necessary to clarify exact details of the swamp
- Climate change impacts are likely to affect DOC flux
- Longer growing season and increased precipitation could lead to **increased DOC concentration, discharge, and flux**
- Higher frequency of intense storms furthers **importance of large DOC “pulse” events** and high-resolution data collection




Fig. 10 DOC flux plotted against discharge for each flux observation.

Acknowledgements

I would like to thank my mentors, Abby Whittington and Professor Cédric G. Fichot, for their invaluable guidance and support throughout this project, and the RISE program for this excellent research opportunity.

References

- Boose E, VanScoy M. 2024. Prospect Hill Hydrological Stations at Harvard Forest since 2005. Harvard Forest Data Archive: HF070 (v.35). Environmental Data Initiative: <https://doi.org/10.6073/pasta/8cd773bd8046f42e5d4bb8826fce6f24>.
- Finzi, A. C., Giasson, M., Barker Plotkin, A. A., Aber, J. D., Boose, E. R., Davidson, E. A., ... Foster, D. R. (2020a). Carbon budget of the Harvard Forest Long-Term ecological research site: Pattern, process, and response to global change. Ecological Monographs, 90(4). doi:10.1002/ecm.1423